

Product datasheet for **TA369863**

C17orf66 (HEATR9) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 50-300 Positive control: Human thyroid cancer Predicted cell location: Nucleus
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Full length fusion protein
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	HEAT repeat containing 9
Database Link:	Entrez Gene 256957 Human A2RTY3
Background:	The function of this protein remains unknown.




[View online »](#)

Product images:



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using TA369863 (HEATR9 Antibody) at dilution 1/110 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using TA369863 (HEATR9 Antibody) at dilution 1/110, treated with fusion protein. (Original magnification: $\times 200$)