

Product datasheet for **TA369873S**

C1orf112 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 100-300 Positive control: Human liver cancer Predicted cell location: Cytoplasm
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human C1orf112
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	chromosome 1 open reading frame 112
Database Link:	Entrez Gene 55732 Human Q9NSG2
Background:	The function of this protein remains unknown.
Synonyms:	FLJ10706; FLJ13470; MGC130018; MGC130019; RP1-97P20.1




[View online »](#)

Product images:



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA369873] (C1orf112 Antibody) at dilution 1/110 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA369873] (C1orf112 Antibody) at dilution 1/110, treated with fusion protein. (Original magnification: ×200)