

## Product datasheet for **TA369916**

### CCDC113 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 50-300 Positive control: Human liver cancer Predicted cell location: Cytoplasm or Nucleus
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human CCDC113
Formulation:	pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	coiled-coil domain containing 113
Database Link:	<a href="#">Entrez Gene 29070 Human Q9H0I3</a>
Background:	Component of centriolar satellites contributing to primary cilium formation.
Synonyms:	DKFZp434N1418; FLJ23555; HSPC065




[View online »](#)


## Product images:




Immunohistochemistry of paraffin-embedded Human liver cancer tissue using TA369916 (CCDC113 Antibody) at dilution 1/65 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using TA369916 (CCDC113 Antibody) at dilution 1/65, treated with fusion protein. (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using TA369916 (CCDC113 Antibody) at dilution 1/65 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using TA369916 (CCDC113 Antibody) at dilution 1/65, treated with fusion protein. (Original magnification: ×200)