

Product datasheet for **TA369997S**


CMPK1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 200-1000 WB positive control: HT-29□Hela□HEPG2□A549□Jurkat and 231 cell lysate IHC: 20-100 Positive control: Human tonsil Predicted cell location: Cytoplasm or Nucleus
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Full length fusion protein
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	22 kDa
Gene Name:	cytidine/uridine monophosphate kinase 1
Database Link:	Entrez Gene 51727 Human P30085
Background:	This gene encodes one of the enzymes required for cellular nucleic acid biosynthesis. This enzyme catalyzes the transfer of a phosphate group from ATP to CMP, UMP, or dCMP, to form the corresponding diphosphate nucleotide. Alternate splicing results in both coding and non-coding transcript variants.
Synonyms:	CMK; CMPK; UCK; UMK; UMP-CMPK; UMP/CMPK; UMPK



[View online »](#)

Product images:

Gel: 12%SDS-PAGE

Lysate: 40 μ g

Lane 1-6: HT-29

Hela

HEPG2


A549

Jurkat and 231 cell lysate


Primary antibody: [TA369997] (CMPK1 Antibody)
at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at
1/8000 dilution


Exposure time: 10 seconds




Immunohistochemistry of paraffin-embedded
Human tonsil tissue using [TA369997] (CMPK1
Antibody) at dilution 1/20 (Original magnification:
 $\times 200$)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA369997] (CMPK1 Antibody) at dilution 1/20, treated with fusion protein. (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using [TA369997] (CMPK1 Antibody) at dilution 1/20 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using [TA369997] (CMPK1 Antibody) at dilution 1/20, treated with fusion protein. (Original magnification: $\times 200$)