


Product datasheet for **TA370043S**

CWF19L1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 20–100 Positive control: Human cervical cancer Predicted cell location: Cytoplasm and Nucleus
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human CWF19L1
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at –20°C.
Stability:	1 year
Gene Name:	CWF19-like 1, cell cycle control (S. pombe)
Database Link:	Entrez Gene 55280 Human Q69YN2
Background:	This gene encodes a member of the CWF19 protein family. Mutations in this gene have been associated with autosomal recessive spinocerebellar ataxia–17 and mild cognitive disability. Alternative splicing results in multiple transcript variants.
Synonyms:	FLJ10998




Product images:


Immunohistochemistry of paraffin-embedded Human cervical cancer tissue using [TA370043] (CWF19L1 Antibody) at dilution 1/20 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human cervical cancer tissue using [TA370043] (CWF19L1 Antibody) at dilution 1/20, treated with fusion protein. (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using [TA370043] (CWF19L1 Antibody) at dilution 1/20 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using [TA370043] (CWF19L1 Antibody) at dilution 1/20, treated with fusion protein. (Original magnification: $\times 200$)