

Product datasheet for **TA370067S**


DDA1 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 20-100 Positive control: Human brain Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Full length fusion protein
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	DET1 and DDB1 associated 1
Database Link:	Entrez Gene 79016 Human Q9BW61
Background:	May be involved in ubiquitination and subsequent proteasomal degradation of target proteins. Component of the DDD-E2 complexes which may provide a platform for interaction with CUL4A and WD repeat proteins.
Synonyms:	C19orf58; MGC2594; PCIA-1; PCIA1



[View online »](#)

Product images:

Immunohistochemistry of paraffin-embedded Human brain tissue using [TA370067] (DDA1 Antibody) at dilution 1/20 (Original magnification: x200)



Immunohistochemistry of paraffin-embedded Human brain tissue using [TA370067] (DDA1 Antibody) at dilution 1/20, treated with fusion protein. (Original magnification: x200)