

## Product datasheet for **TA370112S**


### HAS3 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 500-2000 WB positive control: Human cervical cancer tissue and A549 cell lysates IHC: 25-100 Positive control: Human cervical cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human HAS3
Formulation:	pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	63 kDa
Gene Name:	hyaluronan synthase 3
Database Link:	<a href="#">Entrez Gene 3038 Human O00219</a>
Background:	The protein encoded by this gene is involved in the synthesis of the unbranched glycosaminoglycan hyaluronan, or hyaluronic acid, which is a major constituent of the extracellular matrix. This gene is a member of the NODC/HAS gene family. Compared to the proteins encoded by other members of this gene family, this protein appears to be more of a regulator of hyaluronan synthesis. Alternative splicing results in multiple transcript variants.
Synonyms:	HAS3



[View online »](#)

**Product images:**


Gel: 8%SDS-PAGE


Lysate: 40  $\mu$ g

Lane 1-2: Human cervical cancer tissue and A549 cell lysates


Primary antibody: [TA370112] (HAS3 Antibody) at dilution 1/400

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 10 minutes



Immunohistochemistry of paraffin-embedded Human cervical cancer tissue using [TA370112] (HAS3 Antibody) at dilution 1/20 (Original magnification:  $\times$ 200)



Immunohistochemistry of paraffin-embedded Human cervical cancer tissue using [TA370112] (HAS3 Antibody) at dilution 1/20, treated with fusion protein. (Original magnification:  $\times$ 200)