

Product datasheet for **TA370116**


MPPED2 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 30-150 Positive control: Human tonsil Predicted cell location: Cytoplasm and Cell membrane
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Full length fusion protein
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	metallophosphoesterase domain containing 2
Database Link:	Entrez Gene 744 Human Q15777
Background:	This gene likely encodes a metallophosphoesterase. The encoded protein may play a role a brain development. Alternatively spliced transcript variants have been described.
Synonyms:	239FB; C11orf8; D11S302E; dj873F21.1; dj1024C24.1; FAM1B



[View online »](#)

Product images:

Immunohistochemistry of paraffin-embedded Human tonsil tissue using TA370116 (MPPED2 Antibody) at dilution 1/25 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using TA370116 (MPPED2 Antibody) at dilution 1/25, treated with fusion protein. (Original magnification: ×200)