

Product datasheet for **TA370141S**

MRPL3 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 500-2000 WB positive control: A549, HeLa, HepG2, 293T cell lysates IHC: 50-200 Positive control: Human liver cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human MRPL3
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	39 kDa
Gene Name:	mitochondrial ribosomal protein L3
Database Link:	Entrez Gene 11222 Human P09001




[View online »](#)

Background:


Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 39S subunit protein that belongs to the L3P ribosomal protein family. A pseudogene corresponding to this gene is found on chromosome 13q.

Synonyms:


L3mt; MRL3; MRP-L3; RPML3

Product images:


Gel: 8%SDS-PAGE
 Lysate: 40 µg
 Lane 1-4: A549
 HepG2
 293T cell lysates
 Primary antibody: [TA370141] (MRPL3 Antibody)
 at dilution 1/650
 Secondary antibody: Goat anti rabbit IgG at
 1/5000 dilution
 Exposure time: 30 seconds



Immunohistochemistry of paraffin-embedded
 Human liver cancer tissue using [TA370141]
 (MRPL3 Antibody) at dilution 1/50 (Original
 magnification: x200)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA370141] (MRPL3 Antibody) at dilution 1/50, treated with fusion protein. (Original magnification: ×200)