

Product datasheet for **TA370163S**

DNAJC14 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 50-300 Positive control: Human tonsil Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human DNAJC14
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	Dnaj heat shock protein family (Hsp40) member C14
Database Link:	Entrez Gene 85406 Human Q6Y2X3
Background:	Regulates the export of target proteins, such as DRD1, from the endoplasmic reticulum to the cell surface.
Synonyms:	DNAJ; DRIP78; FLJ32792; HDJ-3; HDJ3; LIP6



[View online »](#)

Product images:

Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA370163] (DNAJC14 Antibody) at dilution 1/55 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA370163] (DNAJC14 Antibody) at dilution 1/55, treated with fusion protein. (Original magnification: ×200)