

Product datasheet for **TA370189S**

ZBTB8 (ZBTB8A) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 500-2000 WB positive control: Human thyroid tissue lysate IHC: 50-300 Positive control: Human thyroid cancer Predicted cell location: Nucleus
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human ZBTB8A
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	50 kDa
Gene Name:	zinc finger and BTB domain containing 8A
Database Link:	Entrez Gene 653121 Human Q96BR9
Background:	May be involved in transcriptional regulation.
Synonyms:	BOZ-F1; BOZF1; FLJ90065; MGC17919; ZBTB8




[View online »](#)


Product images:




Gel: 8%SDS-PAGE
Lysate: 40 μ g
Lane: Human thyroid tissue lysate
Primary antibody: [TA370189] (ZBTB8A Antibody) at dilution 1/350
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution
Exposure time: 15 seconds




Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using [TA370189] (ZBTB8A Antibody) at dilution 1/100 (Original magnification: \times 200)



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using [TA370189] (ZBTB8A Antibody) at dilution 1/100, treated with fusion protein. (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using [TA370189] (ZBTB8A Antibody) at dilution 1/100 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using [TA370189] (ZBTB8A Antibody) at dilution 1/100, treated with fusion protein. (Original magnification: $\times 200$)