

## Product datasheet for **TA370327**

### DTYMK Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 500-2000 WB positive control: Human fetal liver tissue, HepG2, K562 and Raji cell lysates
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human DTYMK
Formulation:	pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	24 kDa
Gene Name:	deoxythymidylate kinase
Database Link:	<a href="#">Entrez Gene 1841 Human P23919</a>
Background:	Catalyzes the conversion of dTMP to dTDP.
Synonyms:	CDC8; FLJ44192; MGC198617; PP3731; TMPK; TYMK



[View online »](#)

## Product images:



Gel: 8%SDS-PAGE

Lysate: 40  $\mu$ g

Lane 1-4: Human fetal liver tissue

HepG2

K562 and Raji cell lysates

Primary antibody: TA370327 (DTYMK Antibody) at dilution 1/400

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 5 seconds