

Product datasheet for **TA370375S**


TTLL1 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 40-200 Positive control: Human tonsil Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human TTLL1
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	tubulin tyrosine ligase like 1
Database Link:	Entrez Gene 25809 Human O95922
Background:	Catalytic subunit of the neuronal tubulin polyglutamylase complex. Modifies alpha- and beta-tubulin, generating side chains of glutamate on the gamma-carboxyl groups of specific glutamate residues within the C-terminal tail of alpha- and beta-tubulin (By similarity).
Synonyms:	C22orf7; HS323M22B; KIAA0173; PGs3



[View online »](#)

Product images:


Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA370375] (TTLL1 Antibody) at dilution 1/60 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA370375] (TTLL1 Antibody) at dilution 1/60, treated with fusion protein. (Original magnification: ×200)