

Product datasheet for **TA370377S**

MRPS5 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 500-2000 WB positive control: A549, HepG2, Raji cell, Rat kidney tissue lysates IHC: 50-100 Positive control: Human colorectal cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human MRPS5
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	48 kDa
Gene Name:	mitochondrial ribosomal protein S5
Database Link:	Entrez Gene 64969 Human P82675




[View online »](#)

Background:


Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 28S subunit protein that belongs to the ribosomal protein S5P family. Pseudogenes corresponding to this gene are found on chromosomes 4q, 5q, and 18q.

Synonyms:


MRP-S5; S5mt

Product images:



Gel: 8%SDS-PAGE
 Lysate: 40 µg
 Lane 1-4: A549
 HepG2
 Raji cell
 Rat kidney tissue lysates
 Primary antibody: [TA370377] (MRPS5 Antibody)
 at dilution 1/700
 Secondary antibody: Goat anti rabbit IgG at
 1/5000 dilution
 Exposure time: 1 minute




Immunohistochemistry of paraffin-embedded
 Human colorectal cancer tissue using [TA370377]
 (MRPS5 Antibody) at dilution 1/75 (Original
 magnification: x200)



Immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using [TA370377] (MRPS5 Antibody) at dilution 1/75, treated with fusion protein. (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA370377] (MRPS5 Antibody) at dilution 1/75 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA370377] (MRPS5 Antibody) at dilution 1/75, treated with fusion protein. (Original magnification: $\times 200$)