

## Product datasheet for **TA370426S**

### WDR47 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 30-150 Positive control: Human tonsil Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human WDR47
Formulation:	pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	WD repeat domain 47
Database Link:	<a href="#">Entrez Gene 22911 Human O94967</a>
Background:	The function of this protein remains unknown.
Synonyms:	FLJ90135; KIAA0893




[View online »](#)

## Product images:



Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA370426] (WDR47 Antibody) at dilution 1/35 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA370426] (WDR47 Antibody) at dilution 1/35, treated with fusion protein. (Original magnification: ×200)