

## Product datasheet for **TA370503S**


### ZNF581 Rabbit Polyclonal Antibody

#### Product data:


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 50-300 Positive control: Human tonsil Predicted cell location: Cytoplasm
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human ZNF581
Formulation:	pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	zinc finger protein 581
Database Link:	<a href="#">Entrez Gene 51545 Human</a> <a href="#">Q9P0T4</a>
Background:	ZNF581 contains 4 C2H2-type zinc fingers. ZNF581 may be involved in transcriptional regulation.
Synonyms:	FLJ22550; HSPC189



[View online »](#)

**Product images:**

Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA370503] (ZNF581 Antibody) at dilution 1/50 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA370503] (ZNF581 Antibody) at dilution 1/50, treated with fusion protein. (Original magnification: ×200)