

## Product datasheet for **TA370513**


### **C1orf83 (TCEANC2) Rabbit Polyclonal Antibody**

#### **Product data:**


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 200-400 Positive control: Human liver cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human TCEANC2
Formulation:	pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	transcription elongation factor A N-terminal and central domain containing 2
Database Link:	<a href="#">Entrez Gene 127428 Human Q96MN5</a>
Background:	The function of this protein remains unknown.
Synonyms:	2210010B22Rik; 2210012G02Rik; 6330404A07Rik; AU015593; Tdeanc2



[View online »](#)

**Product images:**

Immunohistochemistry of paraffin-embedded Human liver cancer tissue using TA370513 (TCEANC2 Antibody) at dilution 1/160 (Original magnification:  $\times 200$ )



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using TA370513 (TCEANC2 Antibody) at dilution 1/160, treated with fusion protein. (Original magnification:  $\times 200$ )