


Product datasheet for **TA370547S**

RNF148 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 50-300 Positive control: Human brain Predicted cell location: Cell membrane
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human RNF148
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	ring finger protein 148
Database Link:	Entrez Gene 378925 Human Q8N7C7
Background:	RNF148 is a single-pass membrane protein. It contains 1 PA (protease associated) domain and 1 RING-type zinc finger. The exact function of RNF148 remains unknown.
Synonyms:	MGC35222; OTTHUMP00000211456

[View online »](#)

Product images:


Immunohistochemistry of paraffin-embedded Human brain tissue using [TA370547] (RNF148 Antibody) at dilution 1/50 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human brain tissue using [TA370547] (RNF148 Antibody) at dilution 1/50, treated with fusion protein. (Original magnification: ×200)