

Product datasheet for **TA370567S**

Rhotekin 2 (RTKN2) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 50-300 Positive control: Human colorectal cancer Predicted cell location: Cytoplasm
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human RTKN2
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	rhotekin 2
Database Link:	Entrez Gene 219790 Human Q8IZC4
Background:	May play an important role in lymphopoiesis.
Synonyms:	bA531F24.1; DKFZp686J10120; Em:AC024597.2; FLJ39352; PLEKHK1




[View online »](#)


Product images:




Immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using [TA370567] (RTKN2 Antibody) at dilution 1/75 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using [TA370567] (RTKN2 Antibody) at dilution 1/75, treated with fusion protein. (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA370567] (RTKN2 Antibody) at dilution 1/75 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA370567] (RTKN2 Antibody) at dilution 1/75, treated with fusion protein. (Original magnification: ×200)