

Product datasheet for **TA370609S**

KPC2 (UBAC1) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 50-300 Positive control: Human tonsil Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human UBAC1
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	UBA domain containing 1
Database Link:	Entrez Gene 10422 Human Q9BSL1
Background:	Non-catalytic subunit of the KPC complex that acts as E3 ubiquitin-protein ligase. Required for poly-ubiquitination and proteasome-mediated degradation of CDKN1B during G1 phase of the cell cycle.
Synonyms:	GBDR1; KPC2; RP11-432J22.3; UBADC1




[View online »](#)

Product images:



Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA370609] (UBAC1 Antibody) at dilution 1/50 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA370609] (UBAC1 Antibody) at dilution 1/50, treated with fusion protein. (Original magnification: ×200)