

Product datasheet for **TA370644**


ZNF256 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 100-500 Positive control: Human tonsil Predicted cell location: Cytoplasm
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human ZNF256
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	zinc finger protein 256
Database Link:	Entrez Gene 10172 Human Q9Y2P7
Background:	Transcriptional repressor that plays a role in cell proliferation. Requires TRIM28 for its activity.
Synonyms:	BMZF-3; BMZF3



[View online »](#)

Product images:

Immunohistochemistry of paraffin-embedded Human tonsil tissue using TA370644 (ZNF256 Antibody) at dilution 1/120 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using TA370644 (ZNF256 Antibody) at dilution 1/120, treated with fusion protein. (Original magnification: $\times 200$)