

## Product datasheet for **TA370671S**

### EEF1D Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 1000-5000 WB positive control: A549, A172 and HepG2 cell lysates IHC: 50-300 Positive control: Human thyroid cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human EEF1D
Formulation:	pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	31 kDa
Gene Name:	eukaryotic translation elongation factor 1 delta
Database Link:	<a href="#">Entrez Gene 1936 Human P29692</a>


**Background:** This gene encodes a subunit of the elongation factor-1 complex, which is responsible for the enzymatic delivery of aminoacyl tRNAs to the ribosome. This subunit, delta, functions as guanine nucleotide exchange factor. It is reported that following HIV-1 infection, this subunit interacts with HIV-1 Tat. This interaction results in repression of translation of host cell proteins and enhanced translation of viral proteins. Several alternatively spliced transcript variants encoding multiple isoforms have been found for this gene. Related pseudogenes have been defined on chromosomes 1, 6, 7, 9, 11, 13, 17, 19.




[View online »](#)

Synonyms: EF-1-delta; EF-1D; EF1D; FLJ20897; FP1047


### Product images:




Gel: 8%SDS-PAGE  
Lysate: 40  $\mu$ g  
Lane 1-3: A549  
A172 and HepG2 cell lysates  
Primary antibody: [TA370671] (EEF1D Antibody)  
at dilution 1/950  
Secondary antibody: Goat anti rabbit IgG at  
1/5000 dilution  
Exposure time: 3 seconds




Immunohistochemistry of paraffin-embedded  
Human thyroid cancer tissue using [TA370671]  
(EEF1D Antibody) at dilution 1/65 (Original  
magnification:  $\times$ 200)



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using [TA370671] (EEF1D Antibody) at dilution 1/65, treated with fusion protein. (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using [TA370671] (EEF1D Antibody) at dilution 1/65 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using [TA370671] (EEF1D Antibody) at dilution 1/65, treated with fusion protein. (Original magnification: ×200)