

Product datasheet for **TA370750**

VILIP3 (HPCAL1) Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 1000-5000 WB positive control: Rat brain tissue, Hela, HepG2 and A172 cell lysates IHC: 50-200 Positive control: Human liver cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human HPCAL1
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	22 kDa
Gene Name:	hippocalcin like 1
Database Link:	Entrez Gene 3241 Human P37235
Background:	The protein encoded by this gene is a member of neuron-specific calcium-binding proteins family found in the retina and brain. It is highly similar to human hippocalcin protein and nearly identical to the rat and mouse hippocalcin like-1 proteins. It may be involved in the calcium-dependent regulation of rhodopsin phosphorylation and may be of relevance for neuronal signalling in the central nervous system. Several alternatively spliced transcript variants encoding the same protein have been found for this gene.




[View online »](#)

Synonyms: BDR1; HLP2; VILIP-3


Product images:




Gel: 12%SDS-PAGE
 Lysate: 40 µg
 Lane 1-4: Rat brain tissue
 HeLa
 HepG2 and A172 cell lysates
 Primary antibody: TA370750 (HPCAL1 Antibody)
 at dilution 1/1000
 Secondary antibody: Goat anti rabbit IgG at
 1/5000 dilution
 Exposure time: 5 seconds




Immunohistochemistry of paraffin-embedded
 Human liver cancer tissue using TA370750
 (HPCAL1 Antibody) at dilution 1/70 (Original
 magnification: ×200)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using TA370750 (HPCAL1 Antibody) at dilution 1/70, treated with fusion protein. (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using TA370750 (HPCAL1 Antibody) at dilution 1/70 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using TA370750 (HPCAL1 Antibody) at dilution 1/70, treated with fusion protein. (Original magnification: ×200)