

Product datasheet for TA370751

CSRP2 Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: IHC, WB

Recommended Dilution: WB: 1000-5000

WB positive control: HepG2 and A549 cell lysates

IHC: 50-200

Positive control: Human tonsil Predicted cell location: Cytoplasm

Reactivity: Human, Mouse, Rat

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

Immunogen: Fusion protein of human CSRP2

Formulation: pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Concentration: lot specific

Purification: Antigen affinity purification

Conjugation: Unconjugated Storage: Store at -20°C.

Stability: 1 year Predicted Protein Size: 21 kDa

Gene Name: cysteine and glycine rich protein 2

Database Link: Entrez Gene 1466 Human

Q16527

Background: CSRP2 is a member of the CSRP family of genes, encoding a group of LIM domain proteins,

which may be involved in regulatory processes important for development and cellular differentiation. CRP2 contains two copies of the cysteine-rich amino acid sequence motif (LIM) with putative zinc-binding activity, and may be involved in regulating ordered cell growth. Other genes in the family include CSRP1 and CSRP3. Alternative splicing results in

multiple transcript variants.



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200


CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com




Synonyms: CRP2; LMO5; SMLIM

Product images:




Gel: 12%SDS-PAGE Lysate: 40 µg Lane 1-2: HepG2 and A549 cell lysates Primary antibody: TA370751 (CSRP2 Antibody) at dilution 1/900 Secondary antibody: Goat anti rabbit IgG at 1/5000 dilution Exposure time: 10 seconds




Immunohistochemistry of paraffin-embedded Human tonsil tissue using TA370751 (CSRP2 Antibody) at dilution 1/60 (Original magnification: ×200)






Immunohistochemistry of paraffin-embedded Human tonsil tissue using TA370751 (CSRP2 Antibody) at dilution 1/60, treated with fusion protein. (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using TA370751 (CSRP2 Antibody) at dilution 1/60 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using TA370751 (CSRP2 Antibody) at dilution 1/60, treated with fusion protein. (Original magnification: ×200)