

Product datasheet for **TA370873S**

NHLRC2 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 150-300 Positive control: Human colorectal cancer Predicted cell location: Cytoplasm
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human NHLRC2
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	NHL repeat containing 2
Database Link:	Entrez Gene 374354 Human Q8NBF2
Background:	Required for normal embryonic development.
Synonyms:	1200003G01Rik; DKFZp779F115; FLJ20147; FLJ25621; FLJ33312; MGC45492




[View online »](#)

Product images:



Immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using [TA370873] (NHLRC2 Antibody) at dilution 1/120 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using [TA370873] (NHLRC2 Antibody) at dilution 1/120, treated with fusion protein. (Original magnification: $\times 200$)