

Product datasheet for **TA371018S**


NOP58 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 500-2000 WB positive control: 293T, HepG2, Hela, 231, K562 cell lysates IHC: 50-200 Positive control: Human gastric cancer Predicted cell location: Nucleus
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human NOP58
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	59 kDa
Gene Name:	NOP58 ribonucleoprotein
Database Link:	Entrez Gene 51602 Human Q9Y2X3
Background:	The protein encoded by this gene is a core component of box C/D small nucleolar ribonucleoproteins. Some box C/D small nucleolar RNAs (snoRNAs), such as U3, U8, and U14, are dependent upon the encoded protein for their synthesis. This protein is SUMOylated, which is necessary for high affinity binding to snoRNAs.
Synonyms:	HSPC120; NOL5; NOP5; NOP5/NOP58


[View online »](#)


Product images:



Gel: 8%SDS-PAGE
 Lysate: 40 μ g
 Lane 1-5: 293T
 HepG2
 Hela
 231
 K562 cell lysates
 Primary antibody: [TA371018] (NOP58 Antibody)
 at dilution 1/600
 Secondary antibody: Goat anti rabbit IgG at
 1/5000 dilution
 Exposure time: 1 second



Immunohistochemistry of paraffin-embedded
 Human gastric cancer tissue using [TA371018]
 (NOP58 Antibody) at dilution 1/85 (Original
 magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human gastric cancer tissue using [TA371018] (NOP58 Antibody) at dilution 1/85, treated with fusion protein. (Original magnification: ×200)