

## Product datasheet for **TA371048S**


### ZNF184 Rabbit Polyclonal Antibody

#### Product data:

|                         |                                                                                                                                                                                                                                                                                                                                                    |
|-------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:           | Primary Antibodies                                                                                                                                                                                                                                                                                                                                 |
| Applications:           | IHC, WB                                                                                                                                                                                                                                                                                                                                            |
| Recommended Dilution:   | WB: 1000-5000<br>WB positive control: RAMOS, Raji, A375 cell, Rat liver tissue, LO2, 293T, Jurkat cell lysates<br>IHC: 100-300<br>Positive control: Human cervical cancer<br>Predicted cell location: Nucleus                                                                                                                                      |
| Reactivity:             | Human, Rat                                                                                                                                                                                                                                                                                                                                         |
| Host:                   | Rabbit                                                                                                                                                                                                                                                                                                                                             |
| Isotype:                | IgG                                                                                                                                                                                                                                                                                                                                                |
| Clonality:              | Polyclonal                                                                                                                                                                                                                                                                                                                                         |
| Immunogen:              | Fusion protein of human ZNF184                                                                                                                                                                                                                                                                                                                     |
| Formulation:            | pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol                                                                                                                                                                                                                                                                                                   |
| Purification:           | Antigen affinity purification                                                                                                                                                                                                                                                                                                                      |
| Conjugation:            | Unconjugated                                                                                                                                                                                                                                                                                                                                       |
| Storage:                | Store at -20°C.                                                                                                                                                                                                                                                                                                                                    |
| Stability:              | 1 year                                                                                                                                                                                                                                                                                                                                             |
| Predicted Protein Size: | 86 kDa                                                                                                                                                                                                                                                                                                                                             |
| Gene Name:              | zinc finger protein 184                                                                                                                                                                                                                                                                                                                            |
| Database Link:          | <a href="#">Entrez Gene 7738 Human Q99676</a>                                                                                                                                                                                                                                                                                                      |
| Background:             | The protein encoded by this gene is predicted to be a Kruppel C2H2-type zinc-finger protein family member. Sequence analysis predicts that the protein contains two Kruppel associated box (KRAB) boxes in the N-terminus and highly conserved zinc finger motifs at the C-terminus. Alternative splicing results in multiple transcript variants. |
| Synonyms:               | OTTHUMP00000016181                                                                                                                                                                                                                                                                                                                                 |



[View online »](#)

**Product images:**


Gel: 6%SDS-PAGE

Lysate: 40 µg

Lane 1-7: RAMOS

Raji

A375 cell

Rat liver tissue

LO2


293T

Jurkat cell lysates


Primary antibody: [TA371048] (ZNF184 Antibody)  
at dilution 1/1000

Secondary antibody: Goat anti rabbit IgG at  
1/5000 dilution

Exposure time: 2 minutes



Immunohistochemistry of paraffin-embedded  
Human cervical cancer tissue using [TA371048]  
(ZNF184 Antibody) at dilution 1/105 (Original  
magnification: ×200)



Immunohistochemistry of paraffin-embedded Human cervical cancer tissue using [TA371048] (ZNF184 Antibody) at dilution 1/105, treated with fusion protein. (Original magnification: ×200)