

## Product datasheet for **TA371287S**


### SIGLEC6 Rabbit Polyclonal Antibody

#### Product data:


Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 500-2000 WB positive control: Human plasma solution IHC: 50-300 Positive control: Human liver cancer Predicted cell location: Secreted and Cell membrane
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide of human SIGLEC6
Formulation:	pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	50 kDa
Gene Name:	sialic acid binding Ig like lectin 6
Database Link:	<a href="#">Entrez Gene 946 Human O43699</a>
Background:	This gene encodes a member of the SIGLEC (sialic acid binding immunoglobulin-like lectin) family of proteins. The encoded transmembrane receptor binds sialyl-TN glycans and leptin. Placental expression of the encoded protein is upregulated in preeclampsia.
Synonyms:	CD33L; CD33L1; CD327; CDw327; OB-BP1; OBBP1; SIGLEC-6




[View online »](#)

**Product images:**

Gel: 8%SDS-PAGE  
Lysate: 40  $\mu$ g  
Lane: Human plasma solution  
Primary antibody: [TA371287] (SIGLEC6 Antibody) at dilution 1/200  
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution  
Exposure time: 40 seconds



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA371287] (SIGLEC6 Antibody) at dilution 1/75 (Original magnification:  $\times$ 200)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA371287] (SIGLEC6 Antibody) at dilution 1/75, treated with synthetic peptide. (Original magnification:  $\times 200$ )