

Product datasheet for **TA371438S**

NAP1 (NAA25) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 50-200 Positive control: Human thyroid cancer Predicted cell location: Cytoplasm
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide of human NAA25
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	N(alpha)-acetyltransferase 25, NatB auxiliary subunit
Database Link:	Entrez Gene 80018 Human Q14CX7
Background:	This gene encodes the auxiliary subunit of the heteromeric N-terminal acetyltransferase B complex. This complex acetylates methionine residues that are followed by acidic or asparagine residues.
Synonyms:	4833422K13Rik; A630046C11; A1450868; C330023M02Rik; LRRGT00164; Mdm20; RGD1305685




[View online »](#)

Product images:



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using [TA371438] (NAA25 Antibody) at dilution 1/40 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using [TA371438] (NAA25 Antibody) at dilution 1/40, treated with synthetic peptide. (Original magnification: ×200)