

Product datasheet for **TA371781S**

DOK6 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 20-100 Positive control: Human ovarian cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide of human DOK6
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	docking protein 6
Database Link:	Entrez Gene 220164 Human Q6PKX4
Background:	DOK6 is a member of the DOK (see DOK1; MIM 602919) family of intracellular adaptors that play a role in the RET (MIM 164761) signaling cascade (Crowder et al., 2004 [PubMed 15286081]).
Synonyms:	DOK5L; HsT3226; MGC20785




[View online »](#)

Product images:



Immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using [TA371781] (DOK6 Antibody) at dilution 1/20 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using [TA371781] (DOK6 Antibody) at dilution 1/20, treated with synthetic peptide. (Original magnification: $\times 200$)