

Product datasheet for **TA371968S**

PLAG1 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 500-2000 WB positive control: Hepg2□A431□Hela and A549 cell lysates IHC: 25-100 Positive control: Human tonsil Predicted cell location: Nucleus
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide of human PLAG1
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	56 kDa
Gene Name:	PLAG1 zinc finger
Database Link:	Entrez Gene 5324 Human Q6DJT9
Background:	Pleomorphic adenoma gene 1 encodes a zinc finger protein with 2 putative nuclear localization signals. PLAG1, which is developmentally regulated, has been shown to be consistently rearranged in pleomorphic adenomas of the salivary glands. PLAG1 is activated by the reciprocal chromosomal translocations involving 8q12 in a subset of salivary gland pleomorphic adenomas. Three transcript variants encoding two different isoforms have been found for this gene.



[View online »](#)

Synonyms: PSA; SGPA

Product images:



Gel: 8%SDS-PAGE

Lysate: 40 μ g

Lane 1-4: Hepg2


A431

Hela and A549 cell lysates


Primary antibody: [TA371968] (PLAG1 Antibody)
at dilution 1/250

Secondary antibody: Goat anti rabbit IgG at
1/8000 dilution

Exposure time: 3 seconds



Immunohistochemistry of paraffin-embedded
Human tonsil tissue using [TA371968] (PLAG1
Antibody) at dilution 1/25 (Original magnification:
 \times 200)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA371968] (PLAG1 Antibody) at dilution 1/25, treated with synthetic peptide. (Original magnification: $\times 200$)