

Product datasheet for **TA372195S**

B7H3 (CD276) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 200-1000 WB positive control: PC-3□231□LNCAP and 293T cell lysates IHC: 25-100 Positive control: Human liver cancer Predicted cell location: Cytoplasm and Cell membrane
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide of human CD276
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	57 kDa
Gene Name:	CD276 molecule
Database Link:	Entrez Gene 80381 Human Q5ZPR3




[View online »](#)

Background:


The protein encoded by this gene belongs to the immunoglobulin superfamily, and thought to participate in the regulation of T-cell-mediated immune response. Studies show that while the transcript of this gene is ubiquitously expressed in normal tissues and solid tumors, the protein is preferentially expressed only in tumor tissues. Additionally, it was observed that the 3' UTR of this transcript contains a target site for miR29 microRNA, and there is an inverse correlation between the expression of this protein and miR29 levels, suggesting regulation of expression of this gene product by miR29. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

Synonyms:


4Ig-B7-H3; B7-H3; B7H3

Product images:


Gel: 8%SDS-PAGE
 Lysate: 40 µg
 Lane 1-4: PC-3
 231
 LNCAP and 293T cell lysates
 Primary antibody: [TA372195] (CD276 Antibody)
 at dilution 1/350
 Secondary antibody: Goat anti rabbit IgG at
 1/8000 dilution
 Exposure time: 15 seconds



Immunohistochemistry of paraffin-embedded
 Human liver cancer tissue using [TA372195]
 (CD276 Antibody) at dilution 1/20 (Original
 magnification: ×200)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA372195] (CD276 Antibody) at dilution 1/20, treated with synthetic peptide. (Original magnification: $\times 200$)