

Product datasheet for **TA372214S**


CED6 (GULP1) Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 20-100 Positive control: Human brain Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide of human GULP1
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	GULP, engulfment adaptor PTB domain containing 1
Database Link:	Entrez Gene 51454 Human Q9UBP9
Background:	The protein encoded by this gene is an adapter protein necessary for the engulfment of apoptotic cells by phagocytes. Several transcript variants, some protein coding and some thought not to be protein coding, have been found for this gene.
Synonyms:	CED-6; CED6; FLJ31156; GULP; OTTHUMP00000205908; OTTHUMP00000205909; OTTHUMP00000205911



[View online »](#)

Product images:

Immunohistochemistry of paraffin-embedded Human brain tissue using [TA372214] (GULP1 Antibody) at dilution 1/35 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human brain tissue using [TA372214] (GULP1 Antibody) at dilution 1/35, treated with synthetic peptide. (Original magnification: ×200)