

Product datasheet for **TA372316S**

PRPF8 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 500-2000 WB positive control: 293T□K562□Raji□Hela and NIH/3T3 cell lysates IHC: 20-100 Positive control: Human esophagus cancer Predicted cell location: Nucleus or Cytoplasm
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide of human PRPF8
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	274 kDa
Gene Name:	pre-mRNA processing factor 8
Database Link:	Entrez Gene 10594 Human Q6P2Q9
Background:	Pre-mRNA splicing occurs in 2 sequential transesterification steps. The protein encoded by this gene is a component of both U2- and U12-dependent spliceosomes, and found to be essential for the catalytic step II in pre-mRNA splicing process. It contains several WD repeats, which function in protein-protein interactions. This protein has a sequence similarity to yeast Prp8 protein. This gene is a candidate gene for autosomal dominant retinitis pigmentosa.
Synonyms:	HPRP8; p220; PRP8; PRPC8; RP13




[View online »](#)


Product images:




Gel: 6%SDS-PAGE
Lysate: 40 μ g
Lane 1-5: 293T
K562
Raji
Hela and NIH/3T3 cell lysates
Primary antibody: [TA372316] (PRPF8 Antibody)
at dilution 1/600
Secondary antibody: Goat anti rabbit IgG at
1/8000 dilution
Exposure time: 50 seconds




Immunohistochemistry of paraffin-embedded
Human esophagus cancer tissue using
[TA372316] (PRPF8 Antibody) at dilution 1/20
(Original magnification: \times 200)



Immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using [TA372316] (PRPF8 Antibody) at dilution 1/20, treated with synthetic peptide. (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA372316] (PRPF8 Antibody) at dilution 1/20 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using [TA372316] (PRPF8 Antibody) at dilution 1/20, treated with synthetic peptide. (Original magnification: $\times 200$)