

Product datasheet for **TA372358S**

CDX4 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 200-1000 WB positive control: Hela and HEPG2 cell □ Human fetal brain tissue lysates
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide of human CDX4
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	30 kDa
Gene Name:	caudal type homeobox 4
Database Link:	Entrez Gene 1046 Human O14627




[View online »](#)

Background:

Cdx1, Cdx2 and Cdx4 are members of the caudal-type homeobox family of genes, which are homologs of the *Drosophila* caudal gene required for anterior-posterior regional identity. The proteins encoded by these genes are transcription factors which play an important role in development by regulating the expression of Hox genes. Hox genes play a fundamental role in the development of the vertebrate central nervous system, heart, axial skeleton, limbs, gut, urogenital tract and external genitalia. Cdx4 is a major positive regulator of the expression of all Hox family members. Due to its critical role as a regulator, Cdx4 is a direct target of the canonical Wnt pathway. The loss of Cdx4 can result in the development of an expanded hindbrain, while the overexpression of Cdx4 may cause the hindbrain to lose its distinct segmental features and resemble the spinal cord.

Synonyms:

CDX4

Product images:

Gel: 8%SDS-PAGE
Lysate: 40 μ g
Lane 1-3: HeLa and HEPG2 cell
Human fetal brain tissue lysates
Primary antibody: [TA372358] (CDX4 Antibody) at dilution 1/200
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution
Exposure time: 5 minutes