

Product datasheet for **TA372535**


FAM108B1 (ABHD17B) Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 30-150 Positive control: Human cervical cancer Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide of human ABHD17B
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	abhydrolase domain containing 17B
Database Link:	Entrez Gene 51104 Human Q5VST6
Background:	The function of this protein remains unknown.
Synonyms:	C9orf77; CGI-67; FAM108B1; RP11-409O11.2



[View online »](#)

Product images:

Immunohistochemistry of paraffin-embedded Human cervical cancer tissue using TA372535 (ABHD17B Antibody) at dilution 1/35 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human cervical cancer tissue using TA372535 (ABHD17B Antibody) at dilution 1/35, treated with synthetic peptide. (Original magnification: $\times 200$)