

Product datasheet for **TA372679**


GPR177 (WLS) Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 20-100 Positive control: Human tonsil Predicted cell location: Cytoplasm and Cell membrane
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide of human WLS
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	wntless Wnt ligand secretion mediator
Database Link:	Entrez Gene 79971 Human Q5T9L3
Background:	Regulates Wnt proteins sorting and secretion in a feedback regulatory mechanism. This reciprocal interaction plays a key role in the regulation of expression, subcellular location, binding and organelle-specific association of Wnt proteins. Plays also an important role in establishment of the anterior-posterior body axis formation during development (By similarity).
Synonyms:	C1orf139; EVI; GPR177; mig-14; MRP



[View online »](#)

Product images:

Immunohistochemistry of paraffin-embedded Human tonsil tissue using TA372679 (WLS Antibody) at dilution 1/45 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using TA372679 (WLS Antibody) at dilution 1/45, treated with synthetic peptide. (Original magnification: ×200)