

Product datasheet for **TA372715**


FAM171A2 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 20-100 Positive control: Human tonsil Predicted cell location: Cytoplasm
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide of human FAM171A2
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	family with sequence similarity 171 member A2
Database Link:	Entrez Gene 284069 Human A8MVW0
Background:	The function of this protein remains unknown.
Synonyms:	MGC34829



[View online »](#)

Product images:

Immunohistochemistry of paraffin-embedded Human tonsil tissue using TA372715 (FAM171A2 Antibody) at dilution 1/40 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using TA372715 (FAM171A2 Antibody) at dilution 1/40, treated with synthetic peptide. (Original magnification: ×200)