

Product datasheet for **TA372830S**

KANSL1L Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 25-100 Positive control: Human esophagus cancer Predicted cell location: Cytoplasm
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide of human KANSL1L
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Gene Name:	KAT8 regulatory NSL complex subunit 1 like
Database Link:	Entrez Gene 151050 Human A0AUZ9
Background:	The function of this protein remains unknown.
Synonyms:	C2orf67; MSL1v2




[View online »](#)


Product images:




Immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using [TA372830] (KANSL1L Antibody) at dilution 1/30 (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using [TA372830] (KANSL1L Antibody) at dilution 1/30, treated with synthetic peptide. (Original magnification: $\times 200$)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA372830] (KANSL1L Antibody) at dilution 1/30 (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using [TA372830] (KANSL1L Antibody) at dilution 1/30, treated with synthetic peptide. (Original magnification: ×200)