

Product datasheet for **TA372989**

TRAF3IP3 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 200-1000 WB positive control: Rat liver tissue and 231 cell lysates
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fusion protein of human TRAF3IP3
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Predicted Protein Size:	64 kDa
Gene Name:	TRAF3 interacting protein 3
Database Link:	Entrez Gene 80342 Human Q9Y228
Background:	The gene encodes a protein that mediates cell growth by modulating the c-Jun N-terminal kinase signal transduction pathway. The encoded protein may also interact with a large multi-protein assembly containing the phosphatase 2A catalytic subunit. Alternative splicing results in multiple transcript variants.
Synonyms:	DJ434O14.3; FLJ44151; MGC117354; MGC163289; T3JAM



[View online »](#)

Product images:



Gel: 8%SDS-PAGE

Lysate: 40 μ g

Lane 1-2: Rat liver tissue and 231 cell lysates

Primary antibody: TA372989 (TRAF3IP3 Antibody)
at dilution 1/500

Secondary antibody: Goat anti rabbit IgG at
1/5000 dilution

Exposure time: 5 minutes