

Product datasheet for **TA373083L**

Goat Anti-Rabbit IgG (H+L) HRP Secondary Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Goat Anti-Rabbit IgG (H+L) HRP Secondary Antibody
Applications:	ELISA, ICC, IHC, WB
Recommended Dilution:	EIA and Western blots: 1:5,000-1:100,000. Immunohistochemistry: 1:500-1:5,000. Optimal working dilutions must be determined by end user.
Reactivity:	Rabbit
Host:	Goat
Immunogen:	Rabbit IgG (H+L)
Isotype:	IgG
Formulation:	1mg/ml in PBS with 0.02% thimerosal, 50% glycerol, pH7.4
Purification:	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Conjugation:	HRP
Storage:	Aliquot and store at -20°C for long term. Avoid freeze-thaw cycles. Prepare fresh working dilution every times.



[View online »](#)