

## Product datasheet for **TA373106**

### A4GNT Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB,1:500 - 1:2000
Reactivity:	Human
Modifications:	Unmodified
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-340 of human A4GNT (NP_057245.1).
Formulation:	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid freeze / thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	39kDa
Gene Name:	alpha-1,4-N-acetylglucosaminyltransferase
Database Link:	<a href="#">Entrez Gene 51146 Human Q9UNA3</a>
Background:	This gene encodes a protein from the glycosyltransferase 32 family. The enzyme catalyzes the transfer of N-acetylglucosamine (GlcNAc) to core 2 branched O-glycans. It forms a unique glycan, GlcNAc $\alpha$ 1 $\rightarrow$ 4Gal $\beta$ 1 $\rightarrow$ 6R and is largely associated with the Golgi apparatus membrane.
Synonyms:	alpha-1,4-N-acetylglucosaminyltransferase; alpha4GnT; MGC149493



[View online »](#)

## Product images:



Western blot analysis of extracts of various cell lines, using A4GNT antibody (TA373106) at 1:1000 dilution. | Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. | Lysates/proteins: 25ug per lane. | Blocking buffer: 3% nonfat dry milk in TBST. | Detection: ECL Basic Kit . | Exposure time: 90s.