

Product datasheet for **TA373379S**

AIFM3 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB,1:500 - 1:2000 IHC,1:50 - 1:200
Reactivity:	Human, Mouse, Rat
Modifications:	Unmodified
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human AIFM3 (NP_001018070.1).
Formulation:	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid freeze / thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	65kDa/66kDa
Gene Name:	apoptosis inducing factor, mitochondria associated 3
Database Link:	Entrez Gene 150209 Human Q96NN9
Background:	Induces apoptosis through a caspase dependent pathway. Reduces mitochondrial membrane potential.
Synonyms:	AIFL; FLJ30473; FLJ45137; OTTHUMP00000198535; OTTHUMP00000198536



[View online »](#)

Product images:



Western blot analysis of extracts of various cell lines, using AIFM3 antibody ([TA373379]) at 1:1000 dilution. | Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. | Lysates/proteins: 25ug per lane. | Blocking buffer: 3% nonfat dry milk in TBST. | Detection: ECL Basic Kit . | Exposure time: 90s.