


## Product datasheet for **TA373392S**

### AK7 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, ICC/IF, WB
Recommended Dilution:	WB, 1:500 - 1:2000 IF/ICC, 1:50 - 1:200 ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Reactivity:	Human, Mouse, Rat
Modifications:	Unmodified
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid freeze / thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	83kDa
Gene Name:	adenylate kinase 7
Database Link:	<a href="#">Entrez Gene 122481 Human</a> <a href="#">Q96M32</a>
Background:	This gene encodes a member of the adenylate kinase family of enzymes. The encoded enzyme is a phosphotransferase that catalyzes the reversible phosphorylation of adenine nucleotides. This enzyme plays a role in energy homeostasis of the cell. Alternative splicing results in multiple transcript variants. Mutations in the mouse gene are associated with primary ciliary dyskinesia.
Synonyms:	FLJ32864

[View online »](#)

**Product images:**

Western blot analysis of various lysates using AK7 Rabbit pAb ([TA373392]) at 1:1000 dilution.  
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Basic Kit (RM00020).  
Exposure time: 3s.