

## Product datasheet for **TA373520**

### ANKRD2 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB,1:1000 - 1:2000
Reactivity:	Mouse, Rat
Modifications:	Unmodified
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 30-140 of human ANKRD2 (NP_001123453.1).
Formulation:	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid freeze / thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	36kDa/39kDa
Gene Name:	ankyrin repeat domain 2
Database Link:	<a href="#">Q9GZV1</a>
Background:	This gene encodes a protein that belongs to the muscle ankyrin repeat protein (MARP) family. A similar gene in rodents is a component of a muscle stress response pathway and plays a role in the stretch-response associated with slow muscle function. Alternative splicing results in multiple transcript variants encoding different isoforms.
Synonyms:	ARPP; hArpp; MGC104314



[View online »](#)

## Product images:



Western blot analysis of extracts of various cell lines, using ANKRD2 antibody (TA373520) at 1:1000 dilution. | Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. | Lysates/proteins: 25ug per lane. | Blocking buffer: 3% nonfat dry milk in TBST. | Detection: ECL Basic Kit . | Exposure time: 5s.