

## **Product datasheet for TA373969**

## **Product data:**

**Product Type:** Primary Antibodies

**Applications:** ICC/IF, IHC, WB

**BCL2 Rabbit Polyclonal Antibody** 

Recommended Dilution: WB,1:500 - 1:2000

IHC,1:50 - 1:200 IF,1:50 - 1:200

Reactivity: Human, Mouse, Rat

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Bcl-2

(NP\_000624.2).

Formulation: Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3

**Concentration:** lot specific

**Purification:** Affinity purification

**Conjugation:** Unconjugated

Storage: Store at -20°C. Avoid freeze / thaw cycles.

**Stability:** Shelf life: one year from despatch.

Predicted Protein Size: 22kDa/26kDa

Gene Name: B-cell CLL/lymphoma 2

Database Link: Entrez Gene 596 Human

P10415

**Background:** This gene encodes an integral outer mitochondrial membrane protein that blocks the

apoptotic death of some cells such as lymphocytes. Constitutive expression of BCL2, such as in the case of translocation of BCL2 to Ig heavy chain locus, is thought to be the cause of

follicular lymphoma. Alternative splicing results in multiple transcript variants.

Synonyms: Bcl-2; OTTHUMP00000163680




**OriGene Technologies, Inc.** 9620 Medical Center Drive, Ste 200


CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



## **Product images:**






HCT116


IL-8

Western blot analysis of extracts of various cells, using Bcl-2 antibody (TA373969) at 1:1000 dilution. | Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. | Lysates/proteins: 25ug per lane. | Blocking buffer: 3% nonfat dry milk in TBST. | Detection: ECL Basic Kit . | Exposure time: 60s.

Western blot analysis of extracts of HCT116 cells, using Bcl-2 antibody (TA373969) at 1:1000 dilution.HCT116 cells treated with different concentrations of IL-8.|Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.|Lysates/proteins: 25ug per lane.|Blocking buffer: 3% BSA.|Detection: ECL Basic Kit.|Exposure time: 30s.







Immunofluorescence analysis of C6 cells using Bcl-2 antibody (TA373969) at dilution of 1:100. Blue: DAPI for nuclear staining.

Immunofluorescence analysis of L929 cells using Bcl-2 antibody (TA373969) at dilution of 1:100. Blue: DAPI for nuclear staining.