

Product datasheet for **TA374081**

BRI3 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB,1:500 - 1:2000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant protein of human BRI3.
Formulation:	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid freeze / thaw cycles.
Stability:	Shelf life: one year from despatch.
Gene Name:	brain protein I3
Database Link:	Entrez Gene 25798 Human O95415
Background:	Participates in tumor necrosis factor-alpha (TNF-induced cell death. May be a target of Wnt/beta-catenin signaling in the liver.
Synonyms:	I3; ITM2C; pRGR2



[View online »](#)