

## Product datasheet for **TA375000S**

### COQ10B Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	WB, 1:500 - 1:1000 ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Reactivity:	Mouse, Rat
Modifications:	Unmodified
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid freeze / thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	27kDa
Gene Name:	coenzyme Q10B
Database Link:	<a href="#">Q9H8M1</a>
Background:	Predicted to enable ubiquinone binding activity. Predicted to be involved in cellular respiration and ubiquinone biosynthetic process. Predicted to be located in mitochondrial inner membrane. Predicted to be active in mitochondrion.
Synonyms:	FLJ13448



[View online »](#)

## Product images:



Western blot analysis of various lysates using COQ10B Rabbit pAb (TA375000S) at 1:1000 dilution.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit.

Exposure time: 10s.