

Product datasheet for **TA376106**

FAM32A Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB,1:500 - 1:2000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant protein of human FAM32A.
Formulation:	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid freeze / thaw cycles.
Stability:	Shelf life: one year from despatch.
Gene Name:	family with sequence similarity 32 member A
Database Link:	Entrez Gene 26017 Human Q9Y421
Background:	Isoform 1, but not isoform 2 or isoform 3, may induce G2 arrest and apoptosis. May also increase cell sensitivity to apoptotic stimuli.
Synonyms:	DKFZP586O0120



[View online »](#)