

## Product datasheet for **TA376108**

### FAM3B Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB,1:500 - 1:2000
Reactivity:	Human, Mouse
Modifications:	Unmodified
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 30-235 of human FAM3B (NP_478066.3).
Formulation:	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid freeze / thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	20kDa/25kDa/30kDa
Gene Name:	family with sequence similarity 3 member B
Database Link:	<a href="#">Entrez Gene 54097 Human P58499</a>


**Background:** FAM3B, also known as Pancreatic-derived factor (PANDER), is an islet-specific secreted cytokine specifically expressed at high levels in the islets of Langerhans of the endocrine pancreas. FAM3B can induce apoptosis of alpha and beta cells in a dose- and time-dependent manner. Previous studies showed that FAM3B regulates glucose and lipid metabolism through interaction with liver and endocrine pancreas. FAM3B silencing activates both extrinsic and intrinsic apoptotic pathways. In general, silencing FAM3B promoted p53 phosphorylation and induced p53 accumulation by decreasing Mdm2 expression, which resulted in apoptotic cell death.



[View online »](#)

Synonyms: 2-21; C21orf11; C21orf76; D21M16SJHU19e; ORF9; PANDER; PRED44

### Product images:



Western blot analysis of extracts of various cell lines, using FAM3B antibody (TA376108) at 1:1000 dilution. | Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. | Lysates/proteins: 25ug per lane. | Blocking buffer: 3% nonfat dry milk in TBST. | Detection: ECL Basic Kit . | Exposure time: 5s.